

Automate your Reactors

Büchi Reactor Control

- Modular control system
- ATEX compliant
- Single process and multipurpose
- Turnkey solutions



Büchi Reactor Control

Standard Setup Features


Display Indication of:

- Temperature (up to 4 values)
- Pressure
- Speed (RPM)

Control of Multiple Variables:

- Pressure/Vacuum
- Inertization Cycles
- Temperature (via Thermostat)
- Stirrer switch (on/off)

→ PLUS Logging of all Values



Upgrade Options

All actuators and display indication of:


- pH (or other Digital Electrode Signals)
- Reflux Ratio
- Speed Adjustment
- 2nd Temperature Control (e.g. cooler)
- Dosing Control

Touchscreen Panels Available for

- Up to 22" non-ATEX
- Up to 12" ATEX Zone 2



→ GMP FDA CRF 21 Part 11 Compliant logging of all values



Please contact us for an offer
or for more customized solutions